

## **Test Definition: FSCNE**

Sugar Cane (Saccharum officinarum) IgE

**Reporting Title:** Sugar Cane IgE **Performing Location:** Eurofins Viracor

#### **Specimen Requirements:**

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

### **Specimen Minimum Volume:**

0.5 mL

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	28 days	

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
Z3102	Sugar Cane IgE	Alphanumeric		7723-0
Z3103	Class	Alphanumeric		102615-2

LOINC and CPT codes are provided by the performing laboratory.

#### **Supplemental Report:**

No

#### **CPT Code Information:**

86003



# **Test Definition: FSCNE**

Sugar Cane (Saccharum officinarum) IgE

R	efe	ren	ce	Va	ılu	es:

<0.35 kU/L