

## **Test Definition: FCTUR**

Turmeric (Curcuma longa) IgE

Reporting Title: Turmeric IgE
Performing Location: Eurofins Viracor

#### **Specimen Requirements:**

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

### **Specimen Minimum Volume:**

0.5 mL

| Specimen Type | Temperature              | Time     | Special Container |
|---------------|--------------------------|----------|-------------------|
| Serum         | Refrigerated (preferred) | 28 days  |                   |
|               | Frozen                   | 365 days |                   |
|               | Ambient                  | 28 days  |                   |

#### **Result Codes:**

| Result ID | Reporting Name | Туре         | Unit | LOINC®   |
|-----------|----------------|--------------|------|----------|
| Z3104     | Turmeric IgE   | Alphanumeric |      | 30041-8  |
| Z3105     | Class          | Alphanumeric |      | 102969-3 |

LOINC and CPT codes are provided by the performing laboratory.

#### **Supplemental Report:**

No

#### **CPT Code Information:**

86003



# **Test Definition: FCTUR**

Turmeric (Curcuma longa) IgE

|  | R | efe | rer | nce | Va | lues: |
|--|---|-----|-----|-----|----|-------|
|--|---|-----|-----|-----|----|-------|

<0.35 kU/L