

Reporting Title: Poplar White (Populus alba) IgE**Performing Location:** Eurofins Viracor**Specimen Requirements:**

Draw blood in a plain, red-top tube, serum gel tube are acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume:

0.5 mL

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	28 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z3114	Poplar White (Populus alba) IgE	Alphanumeric		21473-4
Z3115	Class	Alphanumeric		102802-6

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

86003

Reference Values:

<0.35 kU/L