

**Reporting Title:** IL-8, Serum**Performing Location:** Eurofins Viracor**Specimen Requirements:**

Specimen Type: Serum

Collection Container/Tube: Red or SST

Submission Container/Tube: Plastic vial

Specimen Volume: 1 mL

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 1 mL of serum frozen in a plastic vial.

**Specimen Minimum Volume:**

1 mL

Specimen Type	Temperature	Time	Special Container
Serum	Frozen	14 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
FIL8S	IL-8, Serum	Alphanumeric		33211-4

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

83520

**Reference Values:**

&lt;12.4 pg/mL