

**Reporting Title:** Phadiatop (Allergy Screen)**Performing Location:** Eurofins Viracor**Specimen Requirements:**

Specimen Type: Serum

Container/Tube: Red or SST

Specimen Volume: 0.5 mL

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

**Specimen Minimum Volume:**

0.5 mL

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	28 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
FPHAS	Phadiatop (Allergy Screen)	Alphanumeric		57036-6

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

86005

---

**Reference Values:**

Negative