

Reporting Title: Pneumococcal Ab Panel (12 Serotype)**Performing Location:** Eurofins Viracor**Specimen Requirements:**

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume:

0.5 mL

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	7 days	
	Frozen	365 days	
	Ambient	7 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z3132	Pneumo Ab Type 1	Alphanumeric		85955-3
Z3133	Pneumo Ab Type 3	Alphanumeric		86081-7
Z3134	Pneumo Ab Type 4	Alphanumeric		86108-8
Z3135	Pneumo Ab Type 8	Alphanumeric		86148-4
Z3136	Pneumo Ab Type 9 (9N)	Alphanumeric		86166-6
Z3137	Pneumo Ab Type 12 (12F)	Alphanumeric		85974-4
Z3138	Pneumo Ab Type 14	Alphanumeric		85992-6
Z3139	Pneumo Ab Type 19 (19F)	Alphanumeric		86021-3
Z3140	Pneumo Ab Type 23 (23F)	Alphanumeric		86061-9
Z3141	Pneumo Ab Type 26 (6B)	Alphanumeric		40905-2
Z3142	Pneumo Ab Type 51 (7F)	Alphanumeric		40911-0
Z3143	Pneumo Ab Type 56 (18C)	Alphanumeric		40913-6

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

86317 x 12

Reference Values:

Pneumococcal Antibody Panel	Reference Range
Pneumo Ab Type 1	>1.3 ug/mL
Pneumo Ab Type 3	>1.3 ug/mL
Pneumo Ab Type 4	>1.3 ug/mL
Pneumo Ab Type 8	>1.3 ug/mL
Pneumo Ab Type 9 (9N)	>1.3 ug/mL
Pneumo Ab Type 12 (12F)	>1.3 ug/mL
Pneumo Ab Type 14	>1.3 ug/mL
Pneumo Ab Type 19 (19F)	>1.3 ug/mL
Pneumo Ab Type 23 (23F)	>1.3 uh/mL
Pneumo Ab Type 26 (6B)	>1.3 ug/mL
Pneumo Ab Type 51 (7F)	>1.3 ug/mL
Pneumo Ab Type 56 (18C)	>1.3 ug/mL