

**Test Definition: FBSLG** 

Basil IgG

Reporting Title: Basil IgG

Performing Location: Eurofins Viracor

## **Specimen Requirements:**

Specimen Type: Serum Container/Tube: Red Top Specimen Volume: 0.5 mL

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5

mL of serum refrigerated in a plastic vial.

### **Specimen Minimum Volume:**

0.5 mL

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	7 days	

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
FBSLG	Basil IgG	Alphanumeric		63124-2

LOINC and CPT codes are provided by the performing laboratory.

# **Supplemental Report:**

No

### **CPT Code Information:**

86001



**Test Definition: FBSLG** 

Basil IgG

### **Reference Values:**

Reference ranges have not been established for food-specific IgG tests.