

Reporting Title: Antidepressant Screen, Quant, U**Performing Location:** Medtox Laboratories, Inc.**Specimen Requirements:**

Collect 3 mL random urine without preservatives. Ship refrigerated in a plastic container.

Specimen Minimum Volume:

0.6 mL

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	14 days	
	Frozen	30 days	
	Ambient	72 hours	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z3350	Antidepressant Screen	Alphanumeric		32136-4
Z3351	Amitriptyline	Alphanumeric		16225-5
Z3352	Nortriptyline	Alphanumeric		61427-1
Z3353	Imipramine	Alphanumeric		61418-0
Z3354	Desipramine	Alphanumeric		61412-3
Z3355	Sertraline	Alphanumeric		32112-5
Z3356	Desmethylsertraline	Alphanumeric		75233-7
Z3357	Fluoxetine	Alphanumeric		19471-2
Z3358	Norfluoxetine	Alphanumeric		19629-5
Z3359	Clomipramine	Alphanumeric		45144-3
Z3360	Desmethyldclomipramine	Alphanumeric		58013-4
Z3361	Doxepin	Alphanumeric		61415-6
Z3362	Desmethyldoxepin	Alphanumeric		58014-2
Z3363	Maprotiline	Alphanumeric		3737-4
Z3364	Protriptyline	Alphanumeric		4001-4

Result ID	Reporting Name	Type	Unit	LOINC®
Z3365	Trimipramine	Alphanumeric		4084-0
Z3366	Paroxetine	Alphanumeric		27109-8
Z3367	Cyclobenzaprine	Alphanumeric		14061-6

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

80333
80337
80369

Reference Values:

Antidepressant screen includes the analysis for: amitriptyline, clomipramine and desmethylclomipramine, cyclobenzaprine, desipramine, doxepin and desmethyldoxepin, fluoxetine and norfluoxetine, imipramine, maprotiline, nortriptyline, paroxetine, protriptyline, sertraline and desmethylertraline, trimipramine.