

**Test Definition: FIL2M** 

Interleukin 2

Reporting Title: Interleukin 2

Performing Location: ARUP Laboratories

# **Specimen Requirements:**

Serum

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 1 mL of serum frozen in a plastic vial.

Separate specimens must be submitted when multiple tests are ordered.

Note: Cytokine levels may demonstrate diurnal variation. For longitudinal comparison, it is recommended that cytokine levels be determined at the same time of day.

### **Specimen Minimum Volume:**

0.4 mL

Specimen Type	Temperature	Time	Special Container
Varies	Frozen	365 days	

### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
FIL2M	Interleukin 2	Alphanumeric		33939-0

LOINC and CPT codes are provided by the performing laboratory.

### Supplemental Report:

No

### **CPT Code Information:**

83520



**Test Definition: FIL2M** 

Interleukin 2

## **Reference Values:**

2.1 pg/mL or less