

Reporting Title: Anti-IgE**Performing Location:** Eurofins Viracor**Specimen Requirements:**

Specimen Type: Serum

Container/Tube: Red Top

Specimen Volume: 1.0 mL

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 1.0 mL of serum ambient in a plastic vial.

Specimen Minimum Volume:

0.5 mL

Specimen Type	Temperature	Time	Special Container
Serum	Ambient (preferred)	7 days	
	Frozen	28 days	
	Refrigerated	7 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
FANTI	Anti-IgE	Alphanumeric		Not Provided

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

83516

Reference Values:

Normal