

**Reporting Title:** Fig IgE**Performing Location:** Eurofins Viracor**Specimen Requirements:**

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

**Specimen Minimum Volume:**

0.5 mL

| Specimen Type | Temperature              | Time     | Special Container |
|---------------|--------------------------|----------|-------------------|
| Serum         | Refrigerated (preferred) | 28 days  |                   |
|               | Frozen                   | 365 days |                   |
|               | Ambient                  | 28 days  |                   |

**Result Codes:**

| Result ID | Reporting Name | Type         | Unit | LOINC®   |
|-----------|----------------|--------------|------|----------|
| Z3870     | Fig IgE        | Alphanumeric |      | 11173-2  |
| Z3871     | CLASS          | Alphanumeric |      | 102738-2 |

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

86003

**Reference Values:**

&lt;0.35 kU/L