

Poultry and Meat Panel IgG

Reporting Title: Poultry and Meat Panel IgG **Performing Location:** Eurofins Viracor

Specimen Requirements:

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 1 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume:

1 mL

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	7 days	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
Z3961	Beef IgG	Alphanumeric		35529-7
Z3962	Chicken IgG	Alphanumeric		35532-1
Z3963	Egg White IgG	Alphanumeric		35535-4
Z3964	Egg Yolk IgG	Alphanumeric		35536-2
Z3965	Lamb-Mutton IgG	Alphanumeric		60380-3
Z3966	Pork IgG	Alphanumeric		35543-8
Z3967	Turkey IgG	Alphanumeric		60411-6

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

86001 x 7



Reference Values:

Reference ranges have not been established for food-specific IgG tests. The clinical utility of food-specific IgG tests has not been established. These tests can be used in special clinical situations to select foods for evaluation by diet elimination and challenge in patients who have food-related complaints. It should be recognized that the presence of food-specific IgG alone cannot be taken as evidence of food allergy and only indicates immunologic sensitization by the food allergen in question. This test should only be ordered by physicians who recognize the limitations of the test.