

Bean Navy/White (Phaseolus vulgaris) IgE

## Reporting Title: Bean Navy/White IgE

Performing Location: Eurofins Viracor

#### **Specimen Requirements:**

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

#### **Specimen Minimum Volume:**

0.5 mL

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	28 days	

#### Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
Z3992	Bean Navy/White IgE	Alphanumeric		6279-4
Z3993	CLASS	Alphanumeric		102173-2

LOINC and CPT codes are provided by the performing laboratory.

#### Supplemental Report:

No

### **CPT Code Information:**

86003



# **Test Definition: FNBE**

Bean Navy/White (Phaseolus vulgaris) IgE

#### **Reference Values:**

<0.35 kU/L