Test Definition: FOPE

Perch Ocean

Reporting Title: Perch Ocean IgE **Performing Location:** Eurofins Viracor

Specimen Requirements:

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume:

0.5 mL

Specimen Type	Temperature	Time	Special Container	
Serum	Refrigerated (preferred)	28 days		
	Frozen	365 days		
	Ambient	28 days		

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
Z3994	Perch Ocean IgE	Alphanumeric		Not Provided
Z3995	CLASS	Alphanumeric		Not Provided

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

86003



Test Definition: FOPE

Perch Ocean

Reference Values:

<0.35 kU/L