

Reporting Title: Aspen (Populus tremuloides) IgE**Performing Location:** Eurofins Viracor**Specimen Requirements:**

Specimen Type: Serum

Container/Tube: Red Top

Specimen Volume: 0.5 mL

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume:

0.5 mL

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	28 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z4012	Aspen (Populus tremuloides) IgE	Alphanumeric		7101-9
Z4013	CLASS	Alphanumeric		102152-6

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

86003

Reference Values:

<0.35 kU/L