

**Reporting Title:** Parrot Australian Feathers IgE**Performing Location:** Eurofins Viracor**Specimen Requirements:**

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

**Specimen Minimum Volume:**

0.5 mL

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	28 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
Z4052	Parrot Australian Feathers IgE	Alphanumeric		6030-1
Z4053	CLASS	Alphanumeric		102195-5

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

86003

---

**Reference Values:**

<0.35 kU/L