

Reporting Title: Donath Landsteiner**Performing Location:** ARUP Laboratories**Specimen Requirements:**

Specimen Type: Serum

Container/Tube: Red

Specimen Volume: 3 mL

Draw blood in a plain red-top tube. Maintain specimen at 37 degrees C until serum is separated from cells. Spin down and send 3 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume:

2 mL

Specimen Type	Temperature	Time	Special Container
Serum Red	Refrigerated (preferred)	7 days	
	Ambient	72 hours	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
FDLS	Donath Landsteiner	Alphanumeric		18287-3

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

86940, 86941

Reference Values:

Negative