

Targeted Benzodiazepine Screen, Random, Urine

Reporting Title: Targeted Benzodiazepine Screen, U

Performing Location: Rochester

Additional Testing Requirements:

In most cases, no additional testing is needed after the qualitative targeted benzodiazepine test is performed if the parent or metabolites found are consistent with the patients prescribed medications. However, if unexpected benzodiazepine parent or metabolites are found, confirmatory testing can be requested at an additional charge.

Specimen Requirements:

Supplies: Sarstedt 5 mL Aliquot Tube (T914)
Collection Container/Tube: Plastic urine container
Submission Container/Tube: Plastic, 5-mL tube

Specimen Volume: 1 mL Collection Instructions:

1. Collect a random urine specimen.

2. No preservative

Specimen Minimum Volume:

0.5 mL

Forms:

If not ordering electronically, complete, print, and send a Therapeutics Test Request (T831) with the specimen.

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	14 days	
	Frozen	28 days	
	Ambient	72 hours	

Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
LPPB	LPPB	List prescribed benzodiazepines	Plain Text	Yes

Targeted Benzodiazepine Screen, Random, Urine

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
LPPB	List prescribed benzodiazepines	Numeric		29305-0
604871	Alprazolam	Alphanumeric	ng/mL	94116-1
	Also used by tests: TABSU			
604867	Alpha-Hydroxyalprazolam	Alphanumeric	ng/mL	19325-0
	Also used by tests: TABSU			
604891	Alpha-Hydroxyalprazolam Glucuronide	Alphanumeric	ng/mL	94115-3
	Also used by tests: TABSU			
604872	Chlordiazepoxide	Alphanumeric	ng/mL	19385-4
	Also used by tests: TABSU			
604889	Clobazam	Alphanumeric	ng/mL	94114-6
	Also used by tests: TABSU			
604890	N-Desmethylclobazam	Alphanumeric	ng/mL	94113-8
	Also used by tests: TABSU			
604873	Clonazepam	Alphanumeric	ng/mL	19399-5
	Also used by tests: TABSU			
604267	7-aminoclonazepam	Alphanumeric	ng/mL	94112-0
	Also used by tests: TABSU			
604874	Diazepam	Alphanumeric	ng/mL	19443-1
	Also used by tests: TABSU			
604880	Nordiazepam	Alphanumeric	ng/mL	19624-6
	Also used by tests: TABSU			
604875	Flunitrazepam	Alphanumeric	ng/mL	19466-2
	Also used by tests: TABSU			
604866	7-aminoflunitrazepam	Alphanumeric	ng/mL	94111-2
	Also used by tests: TABSU			
604876	Flurazepam	Alphanumeric	ng/mL	19474-6
	Also used by tests: TABSU			
604868	2-Hydroxy Ethyl Flurazepam	Alphanumeric	ng/mL	94110-4
	Also used by tests: TABSU			

Targeted Benzodiazepine Screen, Random, Urine

Result ID	Reporting Name	Туре	Unit	LOINC®
604877	Lorazepam	Alphanumeric	ng/mL	19520-6
	Also used by tests: TABSU			
604878	Lorazepam Glucuronide	Alphanumeric	ng/mL	94109-6
	Also used by tests: TABSU			
604879	Midazolam	Alphanumeric	ng/mL	19585-9
	Also used by tests: TABSU			
604869	Alpha-Hydroxy Midazolam	Alphanumeric	ng/mL	94108-8
	Also used by tests: TABSU			
604881	Oxazepam	Alphanumeric	ng/mL	19638-6
	Also used by tests: TABSU			
604882	Oxazepam Glucuronide	Alphanumeric	ng/mL	94107-0
	Also used by tests: TABSU			
604883	Prazepam	Alphanumeric	ng/mL	19678-2
	Also used by tests: TABSU			
604884	Temazepam	Alphanumeric	ng/mL	19698-0
	Also used by tests: TABSU			
604885	Temazepam Glucuronide	Alphanumeric	ng/mL	94106-2
	Also used by tests: TABSU			
604886	Triazolam	Alphanumeric	ng/mL	19714-5
	Also used by tests: TABSU			
604870	Alpha-Hydroxy Triazolam	Alphanumeric	ng/mL	94105-4
	Also used by tests: TABSU			
604887	Zolpidem	Alphanumeric	ng/mL	94104-7
	Also used by tests: TABSU			
604888	Zolpidem Phenyl-4-Carboxylic acid	Alphanumeric	ng/mL	94103-9
	Also used by tests: TABSU			
604949	Benzodiazepine Interpretation	Alphanumeric		69050-3
	Also used by tests: TABSU			

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

Targeted Benzodiazepine Screen, Random, Urine

Components:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
LPPB	List prescribed benzodiazepines			Yes	No
TABSU	Targeted Benzodiazepine Screen, U			Yes	No

CPT Code Information:

80347

G0480 (if appropriate)

Reference Values:

Not Detected

Cutoff concentrations: Alprazolam: 10 ng/mL

Alpha-hydroxyalprazolam: 10 ng/mL

Alpha-hydroxyalprazolam glucuronide: 50 ng/mL

Chlordiazepoxide: 10 ng/mL Clobazam: 10 ng/mL

N-desmethylclobazam: 200 ng/mL

Clonazepam: 10 ng/mL

7-Aminoclonazepam: 10 ng/mL

Diazepam: 10 ng/mL Nordiazepam: 10 ng/mL Flunitrazepam: 10 ng/mL 7-Aminoflunitrazepam: 10 ng/mL

Flurazepam: 10 ng/mL

2-Hydroxy ethyl flurazepam: 10 ng/mL

Lorazepam: 10 ng/mL

Lorazepam glucuronide: 50 ng/mL

Midazolam: 10 ng/mL

Alpha-hydroxy midazolam: 10 ng/mL

Oxazepam: 10 ng/mL

Oxazepam glucuronide: 50 ng/mL

Prazepam: 10 ng/mL Temazepam: 10 ng/mL

Temazepam glucuronide: 50 ng/mL

Triazolam: 10 ng/mL

Alpha-hydroxy triazolam: 10 ng/mL

Zolpidem: 10 ng/mL

Zolpidem phenyl-4-carboxylic acid: 10 ng/mL