

## **Test Definition: A1AFS**

Alpha-1-Antitrypsin Clearance, Feces and Serum

Reporting Title: Alpha-1-Antitrypsin Clearance

Performing Location: Rochester

#### **Ordering Guidance:**

The recommended procedure for protein-losing enteropathy is A1AFS / Alpha-1-Antitrypsin Clearance, Feces and Serum.

### **Shipping Instructions:**

Feces and serum should be shipped together. Specimens shipped separately may delay testing.

### **Specimen Requirements:**

Both feces and serum are required.

Blood must be drawn during the stool collection period.

Specimen Type: Serum

Collection Container/Tube: Red top or serum gel

Submission Container/Tube: Plastic vial

Specimen Volume: 1 mL
Collection Instructions:
1. Centrifuge within 2 hours.
2. Aliquot and ship in plastic vial.

Specimen Type: Feces

Supplies: Stool Containers - 24, 48, 72 Hour Kit (T291)

Container/Tube: Stool container Specimen Volume: Entire collection

Collection Instructions:

1. Collect a 24-hour fecal collection.

2. If no specimen is obtained within 24 hours, extend collection time to 48 to 72 hours. Document time frame.

#### **Specimen Minimum Volume:**

Homogenized feces: 1 mL

Serum: 0.5 mL

#### Forms:

If not ordering electronically, complete, print, and send a Gastroenterology and Hepatology Test Request (T728) with the specimen.

# **Test Definition: A1AFS**

Alpha-1-Antitrypsin Clearance, Feces and Serum

| Specimen Type | Temperature        | Time    | Special Container |
|---------------|--------------------|---------|-------------------|
| Fecal         | Frozen (preferred) | 14 days |                   |
|               | Ambient            | 14 days |                   |
|               | Refrigerated       | 14 days |                   |
| Serum         | Frozen (preferred) | 28 days |                   |
|               | Ambient            | 28 days |                   |
|               | Refrigerated       | 28 days |                   |

### Ask at Order Entry (AOE) Questions:

| Test ID | Question ID | Description           | Туре        | Reportable |
|---------|-------------|-----------------------|-------------|------------|
| A1ATF   | HOURS       | Hours: • 24 • 48 • 72 | Answer List | No         |

### **Result Codes:**

| Result ID | Reporting Name                | Туре    | Unit    | LOINC®  |
|-----------|-------------------------------|---------|---------|---------|
| AATS      | Alpha-1-Antitrypsin, S        | Numeric | mg/dL   | 6771-0  |
| AAT24     | Alpha-1-Antitrypsin, 24 Hr, F | Numeric | mg/dL   | 9407-8  |
| CRCLR     | Clearance                     | Numeric | mL/24 h | 18271-7 |

LOINC and CPT codes are provided by the performing laboratory.

## **Supplemental Report:**

No

### **Components:**

| Test ID | Reporting Name                | CPT Units | CPT Code | Always<br>Performed | Orderable<br>Separately |
|---------|-------------------------------|-----------|----------|---------------------|-------------------------|
| AATS    | Alpha-1-Antitrypsin, S        |           |          | Yes                 | No                      |
| A1ATF   | Alpha-1-Antitrypsin, 24 Hr, F |           |          | Yes                 | No                      |

### **CPT Code Information:**



# **Test Definition: A1AFS**

Alpha-1-Antitrypsin Clearance, Feces and Serum

82103 x 2

#### **Reference Values:**

CLEARANCE:

< or =27 mL/24 hours

FECAL ALPHA-1-ANTRYPSIN CONCENTRATION: < or =54 mg/dL

SERUM ALPHA-1-ANTRYPSIN CONCENTRATION: 100-190 mg/dL