

Reporting Title: F2-Isoprostanes, U**Performing Location:** Rochester**Specimen Requirements:**

Patient Preparation: Patient should not take nonsteroidal anti-inflammatory drugs within the 72 hours or aspirin within the 2 weeks prior to specimen collection.

Supplies: Sarstedt 5 mL Aliquot Tube (T914)

Container/Tube: Plastic urine tube

Specimen Volume: 5 mL

Collection Instructions:

1. Collect a random urine specimen.
2. No preservative.

Specimen Minimum Volume:

1 mL

Forms:

If not ordering electronically, complete, print, and send a Cardiovascular Test Request Form (T724) with the specimen.

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	7 days	
	Frozen	30 days	
	Ambient	7 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
88677	15-F2t-Isoprostane, U	Alphanumeric	ng/mg Cr	90783-2

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

82542

Reflex Tests:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
CRETR	Creatinine, Random, U			Yes	No

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
CRETR	CRETR	Creatinine, Random, U	Numeric	mg/dL	2161-8

Reference Values:

> or =18 years: < or =1.0 ng/mg creatinine

Reference values have not been established for patients who are younger than 18 years of age