

### **Test Definition: MCLNG**

MayoComplete Lung Cancer-Targeted Gene Panel with Rearrangement, Tumor

# **Reporting Title:** MayoComplete Lung Cancer Panel **Performing Location:** Rochester

#### Ordering Guidance:

Multiple oncology (cancer) gene panels are available. For more information see Hematology, Oncology, and Hereditary Test Selection Guide.

#### **Necessary Information:**

A pathology report (final or preliminary), at minimum containing the following information, must accompany specimen for testing to be performed:

- 1. Patient name
- 2. Block number-must be on all blocks, slides, and paperwork (can be handwritten on the paperwork)
- 3. Tissue collection date
- 4. Source of the tissue

#### Specimen Requirements:

This assay requires at least 20% tumor nuclei.

-Preferred amount of tumor area with sufficient percent tumor nuclei: tissue 360 mm(2)

-Minimum amount of tumor area: tissue 72 mm(2)

-These amounts are cumulative over up to 15 unstained slides and must have adequate percent tumor nuclei.

-Tissue fixation: 10% neutral buffered formalin, not decalcified

-For specimen preparation guidance, see Tissue Requirements for Solid Tumor Next-Generation Sequencing. In this document, the sizes are given as 4 mm x 4 mm x 10 slides as preferred: approximate/equivalent to 144 mm(2) and the minimum as 3 mm x 1 mm x 10 slides: approximate/equivalent to 36mm(2)

#### Preferred:

Specimen Type: Tissue block

Collection Instructions: Submit a formalin-fixed, paraffin-embedded tissue block with acceptable amount of tumor tissue.

Acceptable:

Specimen Type: Tissue slides

Slides: 1 Stained and 15 unstained

Collection Instructions: Submit 1 slide stained with hematoxylin and eosin and 15 unstained, nonbaked slides wit 5-micron thick sections of the tumor tissue.

Note: The total amount of required tumor nuclei can be obtained by scraping up to 15 slides from the same block. Additional Information: Unused unstained slides will not be returned.

Specimen Type: Cytology slides (direct smears or ThinPrep)

Slides: 2 to 4 Slides

Collection Instructions: Submit 2 to 4 slides stained and coverslipped with a preferred total of 10,000 nucleated cells, or a minimum of at least 3000 nucleated cells.

Note: Glass coverslips are preferred; plastic coverslips are acceptable but will result in longer turnaround times. Additional Information: Cytology slides will not be returned.



### **Test Definition: MCLNG**

MayoComplete Lung Cancer-Targeted Gene Panel with Rearrangement, Tumor

#### **Specimen Minimum Volume:**

See Specimen Required

#### Forms:

If not ordering electronically, complete, print, and send an Oncology Test Request (T729) with the specimen.

Specimen Type	Temperature	Time	Special Container
Varies	Ambient (preferred)		
	Refrigerated		

#### **Result Codes:**

Result ID	Reporting Name	Type Unit		LOINC®
617833	Result	Alphanumeric		82939-0
617834	Interpretation	Alphanumeric		69047-9
617835	Additional Information	Alphanumeric		48767-8
617836	Specimen	Alphanumeric		31208-2
617837	Tissue ID	Alphanumeric		80398-1
617838	Method	Alphanumeric		85069-3
617839	Disclaimer	Alphanumeric		62364-5
617840	Released By	Alphanumeric		18771-6

LOINC and CPT codes are provided by the performing laboratory.

### Supplemental Report:

No

#### **CPT Code Information:**

88381-Microdissection, manual 81457



## **Test Definition: MCLNG**

MayoComplete Lung Cancer-Targeted Gene Panel with Rearrangement, Tumor

#### **Reflex Tests:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
SLIRV	Slide Review in MG			Yes	No (Bill Only)

#### **Reference Values:**

An interpretive report will be provided.