MAYO CLINIC LABORATORIES BRAF and KIT Mutation Analysis, Next-Generation Sequencing, Tumor

Test Definition: BRFKT

Reporting Title: BRAF/KIT Mutation Analysis, Tumor Performing Location: Rochester

Necessary Information:

A pathology report (final or preliminary), at minimum containing the following information, must accompany specimen for testing to be performed:

- 1. Patient name
- 2. Block number-must be on all blocks, slides, and paperwork (can be handwritten on the paperwork)
- 3. Tissue collection date
- Source of the tissue

Specimen Requirements:

This assay requires at least 20% tumor nuclei.

-Preferred amount of tumor area with sufficient percent tumor nuclei: tissue 216 mm(2)

-Minimum amount of tumor area: tissue 36 mm(2)

-These amounts are cumulative over up to 10 unstained slides and must have adequate percent tumor nuclei

-Tissue fixation: 10% neutral buffered formalin, not decalcified

-For specimen preparation guidance, see Tissue Requirement for Solid Tumor Next-Generation Sequencing. In this document, the sizes are given as 4 mm x 4 mm x 10 slides as preferred: approximate/equivalent to 144 mm(2) and the minimum as 3 mm x 1 mm x 10 slides: approximate/equivalent to 36 mm(2).

Preferred:

Specimen Type: Tissue block

Collection Instructions: Submit a formalin-fixed, paraffin-embedded tissue block with acceptable amount of tumor tissue.

Acceptable:

Specimen Type: Tissue slides

Slides: 1 Stained and 10 unstained

Collection Instructions: Submit 1 slide stained with hematoxylin and eosin and 10 unstained, nonbaked slides with 5-micron thick sections of the tumor tissue.

Note: The total amount of required tumor nuclei can be obtained by scraping up to 10 slides from the same block. Additional Information: Unused unstained slides will not be returned.

Specimen Type: Cytology slides (direct smears or ThinPrep)

Slides: 1 to 3 Slides

Collection Instructions: Submit 1 to 3 slides stained and coverslipped with a preferred total of 5000 nucleated cells, or a minimum of at least 3000 nucleated cells.

Note: Glass coverslips are preferred; plastic coverslips are acceptable but will result in longer turnaround times. Additional Information: Cytology slides will not be returned.

Specimen Minimum Volume:

See Specimen Required

Forms:

Test Definition: BRFKT BRAF and KIT Mutation Analysis, Next-Generation Sequencing, Tumor

If not ordering electronically, complete, print, and send an Oncology Test Request (T729) with the specimen.

Specimen Type	Temperature	Time	Special Container
Varies	Ambient (preferred)		
	Refrigerated		

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
617945	Result	Alphanumeric		82939-0
617946	Interpretation	Alphanumeric		69047-9
617947	Additional Information	Alphanumeric		48767-8
617948	Specimen	Alphanumeric		31208-2
617949	Tissue ID	Alphanumeric		80398-1
617950	Method	Alphanumeric		85069-3
617951	Disclaimer	Alphanumeric		62364-5
617952	Released By	Alphanumeric		18771-6

LOINC and CPT codes are provided by the performing laboratory.

MAYO CLINIC LABORATORIES

Supplemental Report:

No

CPT Code Information:

88381-Microdissection, manual 81406 81479 81479 (if appropriate for government payers)

Reflex Tests:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
SLIRV	Slide Review in MG			Yes	No



Reference Values:

An interpretive report will be provided.