



**Reporting Title:** RPR Screen w/ Reflex to Titer,S

**Performing Location:** Rochester

**Ordering Guidance:**

This test is for monitoring response to therapy in patients treated for syphilis infection. This test should not be used as a primary diagnostic approach for syphilis. To screen for undiagnosed syphilis infection, order SYPH1 / Syphilis IgG with Reflex, Enzyme Immunoassay, Serum

**Specimen Requirements:**

Collection Container/Tube:

Preferred: Serum gel

Acceptable: Red top

Submission Container/Tube: Plastic vial

Specimen Volume: 0.5 mL

Collection Instructions: Centrifuge and aliquot serum into a plastic vial.

**Specimen Minimum Volume:**

0.4 mL

**Forms:**

If not ordering electronically, complete, print, and send Infectious Disease Serology Test Request (T916) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	14 days	
	Frozen	14 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
616863	RPR Screen w/ Reflex to Titer,S	Alphanumeric		20507-0

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

86592

86593-Rapid Plasma Reagin Titer (if appropriate)

**Reflex Tests:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
RPRT2	RPR Titer,S			No	No

**Result Codes for Reflex Tests:**

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
RPRT2	616865	RPR Titer,S	Alphanumeric		31147-2

**Reference Values:**

Negative

Reference values apply to all ages.