

Reporting Title: RPR w/ Reflex to TP-PA, S**Performing Location:** Rochester**Ordering Guidance:**

This assay should only be used following an initial reactive treponemal antibody assay result as part of the reverse syphilis testing algorithm.

Specimen Requirements:

Collection Container/Tube:

Preferred: Serum gel

Acceptable: Red top

Submission Container/Tube: Plastic vial

Specimen Volume: 0.3 mL

Collection Instructions: Centrifuge and aliquot serum into a plastic vial.

Specimen Minimum Volume:

See Specimen Required

Forms:

If not ordering electronically, complete, print, and send an Infectious Disease Serology Test Request (T916) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	14 days	
	Frozen	14 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
616970	RPR w/ Reflex to TP-PA, S	Alphanumeric		20507-0

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

86592
86593- (if appropriate)
86780- (if appropriate)

Reflex Tests:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
RPRT4	RPR, Titer, S			No	No
RTPPA	Syphilis Ab, TP-PA, S			No	Yes (Order TPPA)

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
RPRT4	616971	RPR, Titer, S	Alphanumeric		31147-2
RTPPA	34512	Syphilis Ab, TP-PA, S	Alphanumeric		24312-1

Reference Values:

Negative
Reference values apply to all ages.