MAYO CLINIC
 Test Definition: NIVIIIO

 LABORATORIES
 Nuclear Mitochondrial Gene Panel, Next-Generation Sequencing,

Varies

Reporting Title: Nuclear Mitochondrial Gene Panel Performing Location: Rochester

Ordering Guidance:

The diagnostic workup for a mitochondrial disorder may include testing to demonstrate elevations of the lactate-to-pyruvate ratio and an elevated growth differentiation factor 15 concentration. Consider LAPYP / Lactate Pyruvate Panel, Plasma and GDF15 / Growth Differentiation Factor 15, Plasma.

Customization of this panel and single gene analysis for any gene present on this panel are available. For more information see CGPH / Custom Gene Panel, Hereditary, Next-Generation Sequencing, Varies.

Targeted testing for familial variants (also called site-specific or known variants testing) is available for the genes on this panel. See FMTT / Familial Variant, Targeted Testing, Varies. To obtain more information about this testing option, call 800-533-1710.

Shipping Instructions:

Specimen preferred to arrive within 96 hours of collection.

Specimen Requirements:

Patient Preparation: A previous bone marrow transplant from an allogenic donor will interfere with whole blood or dried blood spot testing. For instructions for testing patients who have received a bone marrow transplant, call 800-533-1710.

Submit only 1 of the following specimens:

Specimen Type: Whole blood Container/Tube: Lavender top (EDTA) or yellow top (ACD) Specimen Volume: 3 mL **Collection Instructions:** 1. Invert several times to mix blood. 2. Send whole specimen in original tube. Do not aliquot. Specimen Stability Information: Ambient (preferred) 4 days/Refrigerated 4 days

Specimen Type: Skin biopsy Supplies: Fibroblast Biopsy Transport Media (T115) Container/Tube: Sterile container with any standard cell culture media (eg, minimal essential media, RPMI 1640). The solution should be supplemented with 1% penicillin and streptomycin. Specimen Volume: 4-mm punch Specimen Stability Information: Refrigerated (preferred)/Ambient Additional Information: A separate culture charge will be assessed under CULFB / Fibroblast Culture for Biochemical or Molecular Testing, Chorionic Villi/Products of Conception/Tissue. An additional 3 to 4 weeks is required to culture fibroblasts before genetic testing can occur.

Specimen Type: Cultured fibroblast Container/Tube: T-25 flask Specimen Volume: 2 Flasks

ABORATORIES Nuclear Mitochondrial Gene Panel, Next-Generation Sequencing, Varies

Collection Instructions: Submit confluent cultured fibroblast cells from a skin biopsy from another laboratory. Cultured cells from a prenatal specimen will not be accepted.

Specimen Stability Information: Ambient (preferred)/Refrigerated (<24 hours)

Additional Information: A separate culture charge will be assessed under CULFB / Fibroblast Culture for Biochemical or Molecular Testing, Chorionic Villi/Products of Conception/Tissue. An additional 3 to 4 weeks is required to culture fibroblasts before genetic testing can occur.

Specimen Minimum Volume:

Blood: 1 mL; Skin biopsy or cultured fibroblasts: See Specimen Required

Forms:

1. New York Clients-Informed consent is required. Document on the request form or electronic order that a copy is on file. The following documents are available:

-Informed Consent for Genetic Testing (T576)

-Informed Consent for Genetic Testing (Spanish) (T826)

2. Molecular Genetics: Biochemical Disorders Patient Information (T527)

3. If not ordering electronically, complete, print, and send a Neurology Specialty Testing Client Test Request (T732) with the specimen.

Specimen Type	Temperature	Time	Special Container
Varies	Varies		

Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
CULFB	CG770	Reason for Referral	Plain Text	No
CULFB	CG899	Specimen	Plain Text	No

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
617091	Test Description	Alphanumeric		62364-5
617092	Specimen	Alphanumeric		31208-2
617093	Source	Alphanumeric		31208-2
617094	Result Summary	Alphanumeric		50397-9

 MAYO CLINIC
 Test Definition: NMITO

 LABORATORIES
 Nuclear Mitochondrial Gene Panel, Next-Generation Sequencing,
 Varies

Result ID	Reporting Name	Туре	Unit	LOINC®
617095	Result	Alphanumeric		82939-0
617096	Interpretation	Alphanumeric		69047-9
618172	Additional Results	Alphanumeric		82939-0
617097	Resources	Alphanumeric		99622-3
617098	Additional Information	Alphanumeric		48767-8
617099	Method	Alphanumeric		85069-3
617100	Genes Analyzed	Alphanumeric		48018-6
617101	Disclaimer	Alphanumeric		62364-5
617102	Released By	Alphanumeric		18771-6

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

81440

Reflex Tests:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
CULFB	Fibroblast Culture for Genetic Test			No	Yes

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Туре	Unit	LOINC®
CULFB	52327	Result Summary	Alphanumeric		50397-9
CULFB	52329	Interpretation	Alphanumeric		69965-2
CULFB	52328	Result	Alphanumeric		82939-0
CULFB	CG770	Reason for Referral	Alphanumeric		42349-1
CULFB	CG899	Specimen	Alphanumeric		31208-2



 MAYOCLINIC
 Test Definition: NVIIIO

 LABORATORIES
 Nuclear Mitochondrial Gene Panel, Next-Generation Sequencing, Varies
 Varies

Test ID	Result ID	Reporting Name	Туре	Unit	LOINC®
CULFB	52331	Source	Alphanumeric		31208-2
CULFB	52332	Method	Alphanumeric		85069-3
CULFB	54625	Additional Information	Alphanumeric		48767-8
CULFB	52333	Released By	Alphanumeric		18771-6

Reference Values:

An interpretive report will be provided