

**Reporting Title:** HBZ ISH, Tech Only**Performing Location:** Rochester**Ordering Guidance:**

This test includes only technical performance of the stain (no pathologist interpretation is performed). If diagnostic consultation by a pathologist is required, order PATHC / Pathology Consultation.

**Shipping Instructions:**

Attach the green pathology address label and the pink Immunostain Technical Only label included in the kit to the outside of the transport container.

**Specimen Requirements:**

Supplies: Immunostain Technical Only Envelope (T693)

Specimen Type: Tissue

Container/Tube: Immunostain Technical Only Envelope

Preferred: 4 Unstained positively charged glass slide (25- x 75- x 1-mm) per test ordered; sections 4-microns thick.

Acceptable: Formalin-fixed, paraffin-embedded (FFPE) tissue block

Specimen Type	Temperature	Time	Special Container
TECHONLY	Ambient (preferred)		
	Refrigerated		

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
618482	HBZ ISH, Tech Only	Alphanumeric		No LOINC Needed

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

88365-TC, Primary

88364-TC, if additional ISH

---

**Reflex Tests:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
ISTOI	ISH Initial, Tech Only			No	No
ISTOA	ISH Additional, Tech Only			No	No