

Test Definition: CD200

CD200 (OX2) Immunostain, Technical Component Only

Reporting Title: CD200 (OX2) IHC, Tech Only

Performing Location: Rochester

Ordering Guidance:

This test includes only technical performance of the stain (no pathologist interpretation is performed). If diagnostic consultation by a pathologist is required, order PATHC / Pathology Consultation.

Shipping Instructions:

Attach the green pathology address label and the pink Immunostain Technical Only label included in the kit to the outside of the transport container.

Specimen Requirements:

Supplies: Immunostain Technical Only Envelope (T693)

Specimen Type: Tissue

Container/Tube: Immunostain Technical Only Envelope

Preferred: 2 unstained positively charged glass slide (25- x 75- x 1-mm) per test ordered; sections 4-microns thick.

Acceptable: Formalin-fixed, paraffin-embedded (FFPE) tissue block

Specimen Type	Temperature	Time	Special Container
TECHONLY	Ambient (preferred)		
	Refrigerated		

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
619518	CD200 (OX2) IHC, Tech Only	Alphanumeric		Bill only; no re

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

88342-TC, primary 88341-TC, if additional IHC



Test Definition: CD200

CD200 (OX2) Immunostain, Technical Component Only

Reflex Tests:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
IHTOI	IHC Initial, Tech Only			No	No
IHTOA	OA IHC Additional, Tech Only			No	No