

## **Test Definition: TAPEN**

Tapentadol and Metabolite, Random, Urine

Reporting Title: Tapentadol and Metabolite, U

Performing Location: Rochester

## **Specimen Requirements:**

Supplies: Sarstedt Aliquot Tube, 5 mL (T914) Collection Container/Tube: Plastic urine container

Submission Container/Tube: Plastic vial

Specimen Volume: 1 mL Collection Instructions:

1. Collect a random urine specimen.

2. No preservative.

## **Specimen Minimum Volume:**

0.5 mL

#### Forms:

If not ordering electronically, complete, print, and send a Therapeutics Test Request (T831) with the specimen.

Specimen Type	Temperature	Time	Special Container	
Urine	Refrigerated (preferred)	14 days		
	Ambient	14 days		
	Frozen	14 days		

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
35916	Tapentadol	Alphanumeric	ng/mL	65807-0
35917	N-desmethyltapentadol	Alphanumeric	ng/mL	65808-8

LOINC and CPT codes are provided by the performing laboratory.

## **Supplemental Report:**

No

#### **CPT Code Information:**

80372

G0480 (if appropriate)



# **Test Definition: TAPEN**

Tapentadol and Metabolite, Random, Urine

#### **Reference Values:**

Cutoff: 25 ng/mL

Note: Tapentadol concentrations will be reported quantitatively and N-desmethyltapentadol will be reported qualitatively (Present or Negative).