

Reporting Title: Oxycodone Screen, U**Performing Location:** Rochester**Ordering Guidance:**

For situations where chain of custody is required, a Chain-of-Custody Kit (T282) is available. For chain-of-custody testing, order OXYSX / Oxycodone Screen, Chain of Custody, Random, Urine.

Additional drug panels and specific requests are available; call 800-533-1710 or 507-266-5700.

Additional Testing Requirements:

If urine creatinine is required or adulteration of the sample is suspected, the following test should also be ordered, ADULT / Adulterants Survey, Random, Urine.

Specimen Requirements:

Supplies: Sarstedt 5 mL Aliquot Tube (T914)

Collection Container/Tube: Plastic urine container

Submission Container/Tube: Plastic, 5-mL tube

Specimen Volume: 2 mL

Collection Instructions:

1. Collect a random urine specimen.
2. No preservative.

Specimen Minimum Volume:

0.5 mL

Forms:

If not ordering electronically, complete, print, and send a Therapeutics Test Request (T831) with the specimen.

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	14 days	
	Frozen	14 days	
	Ambient	72 hours	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
62623	Oxycodone	Alphanumeric	ng/mL	19642-8

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

80307

Reflex Tests:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
OXYCU	Oxycodone w/metabolite Conf, U			No	Yes

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
OXYCU	62616	Oxycodone-by LC-MS/MS	Alphanumeric	ng/mL	16249-5
OXYCU	42010	Noroxycodone-by LC-MS/MS	Alphanumeric	ng/mL	61425-5
OXYCU	35965	Oxymorphone-by LC-MS/MS	Alphanumeric	ng/mL	17395-5
OXYCU	42011	Noroxymorphone-by LC-MS/MS	Alphanumeric	ng/mL	90894-7
OXYCU	36022	Oxycodone Interpretation	Alphanumeric		69050-3

Reference Values:

Negative

Screening cutoff concentration:

Oxycodone: 100 ng/mL