

Reporting Title: Cryptococcus Ag w/Reflex, LFA, CSF
Performing Location: Rochester

Specimen Requirements:

Container/Tube: Sterile vial
Specimen Volume: 1 mL

Specimen Minimum Volume:

0.5 mL

Forms:

If not ordering electronically, complete, print, and send Infectious Disease Serology Test Request (T916) with the specimen.

Specimen Type	Temperature	Time	Special Container
CSF	Refrigerated (preferred)	14 days	
	Frozen	14 days	

Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Type	Reportable
FGENC	Q00M0021	Specimen Source	Plain Text	No

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
62074	Cryptococcus Ag Screen w/Titer, CSF Also used by tests: CLFA	Alphanumeric		29896-8

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

Components:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
CLFA	Cryptococcus Ag Screen w/Titer, CSF			Yes	Yes

CPT Code Information:

87899-Cryptococcus Ag Screen w/Titer, CSF
 87899-Cryptococcus Ag Titer, LFA, CSF (as appropriate)
 87102-Fungal Culture, CSF (as appropriate)

Reflex Tests:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
CLFAT	Cryptococcus Ag Titer, LFA, CSF			No	Yes
FGENC	Fungal Culture, CSF			No	Yes (Order FGEN)
D2F	D2 Fungal Sequencing Identification			No	No (Bill Only)
FUNA	Fungal Ident Panel A			No	No (Bill Only)
FUNB	Fungal Ident Panel B			No	No (Bill Only)
LCCI	Ident Rapid PCR Coccidioides			No	No (Bill Only)
LCHB	Id, Histoplasma/Blastomyces PCR			No	No (Bill Only)
RMALF	Id MALDI-TOF Mass Spec Fungi			No	No (Bill Only)
RMALY	Id MALDI-TOF Mass Spec Yeast			No	No (Bill Only)
LCCA	Id, Candida auris Rapid PCR			No	No (Bill Only)

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
CLFAT	62076	Cryptococcus Ag Titer, LFA, CSF	Alphanumeric		9817-8
FGENC	FGENC	Fungal Culture, CSF	Alphanumeric		In Process

Reference Values:**CRYPTOCOCCUS ANTIGEN SCREEN WITH TITER**

Negative

Reference values apply to all ages.

CRYPTOCOCCUS ANTIGEN TITER, LFA

Negative

Reference values apply to all ages.

FUNGAL CULTURE

Negative

If positive, fungus will be identified.

Reference values apply to all ages.