

**Reporting Title:** Drug Scrn, Prescription/OTC, CoC, S**Performing Location:** Rochester**Specimen Requirements:**

Supplies: Chain-of-Custody Kit (T282)

Container/Tube: Chain-of-Custody Kit containing the specimen containers, seals, and documentation required.

Preferred: 10-mL Red top (serum gel/SST are not acceptable)

Acceptable: 5-mL Red top

Submission Container/Tube: Plastic vial

Specimen Volume: 2.75 mL

Collection Instructions:

1. Centrifuge and aliquot serum into plastic vial within 2 hours of collection.
2. Cap, seal, and submit with the associated documentation to satisfy the legal requirements for chain-of-custody testing.

Additional Information: See Emergency/Overdose.

**Specimen Minimum Volume:**

1.1 mL

**Forms:**

1. Chain of Custody Request is included in the Chain-of-Custody Kit (T282).
2. If not ordering electronically, complete, print, and send a Therapeutics Test Request (T831) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum Red	Refrigerated (preferred)	14 days	
	Frozen	14 days	
	Ambient	3 hours	

**Ask at Order Entry (AOE) Questions:**

Test ID	Question ID	Description	Type	Reportable
DSSX	SUDR	Suspect Drug	Plain Text	No

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
36185	Drugs detected:	Alphanumeric		20785-2
36186	Chain of Custody	Alphanumeric		77202-0

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

80307

**Reflex Tests:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
COCH	Chain of Custody Processing			Yes	No

**Reference Values:**

Drugs detected are presumptive. Additional testing may be required to confirm the presence of any drugs detected.