

Reporting Title: Pain Clinic Survey 10, CoC, U**Performing Location:** Rochester**Specimen Requirements:**

Container/Tube: Chain-of-Custody Kit (T282) containing the specimen containers, seals, and documentation required.

Specimen Volume: 30 mL

Collection Instructions: Collect a random specimen in the container provided, seal, and submit with the associated documentation to satisfy the legal requirements for chain-of-custody testing.

Additional Information: Submitting less than 30 mL will compromise our ability to perform all necessary testing.

Specimen Minimum Volume:

20 mL

Forms:

1. Chain of Custody Request is included in the Chain-of-Custody Kit (T282).
2. If not ordering electronically, complete, print, and send a Therapeutics Test Request (T831) with the specimen.

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	7 days	
	Frozen	14 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
36253	Amphetamines Also used by tests: CDA5X, CDAUX, CDA7X	Alphanumeric	ng/mL	43983-6
36258	Barbiturates Also used by tests: CDAUX, CDA7X	Alphanumeric	ng/mL	70155-7
36259	Benzodiazepines Also used by tests: CDAUX, CDA7X	Alphanumeric	ng/mL	14316-4
36254	Cocaine Also used by tests: CDA5X, CDA7X, CDAUX	Alphanumeric	ng/mL	43984-4

Result ID	Reporting Name	Type	Unit	LOINC®
36260	Methadone metabolite Also used by tests: CDA7X	Alphanumeric	ng/mL	41858-2
36255	Opiates Also used by tests: CDA5X, CDA7X, CDAUX	Alphanumeric	ng/mL	70151-6
36256	Phencyclidine Also used by tests: CDA5X, CDA7X, CDAUX	Alphanumeric	ng/mL	14310-7
36257	Tetrahydrocannabinol Also used by tests: CDA7X, CDA5X, CDAUX	Alphanumeric	ng/mL	14312-3
36261	Chain of Custody Also used by tests: CDAUX, CDA5X, CDA7X	Alphanumeric		77202-0
61727	Oxycodone Also used by tests: OXYSX	Alphanumeric	ng/mL	19642-8
36027	Chain of Custody Also used by tests: OXYSX	Alphanumeric		77202-0

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

Components:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
PN10X	Pain Clinic Survey 10, CoC			Yes	No
OXYSX	Oxycodone Screen, CoC, U			Yes	Yes

CPT Code Information:

80307

Reflex Tests:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
COCH	Chain of Custody Processing			Yes	No
ADLTX	Adulterants Survey, CoC, U			Yes	Yes
AMPHX	Amphetamines Confirmation, CoC, U			No	Yes
COKEX	Cocaine and metabolite Conf, CoC, U			No	Yes
BARBX	Barbiturates Confirmation, CoC, U			No	Yes
MTDNX	Methadone Confirmation, CoC, U			No	Yes
PCPX	Phencyclidine Confirmation, CoC, U			No	Yes
THCX	Carboxy-THC Confirmation, CoC, U			No	Yes
OPATX	Opiate Confirmation, CoC, U			No	Yes
OXYCX	Oxycodone w/metabolite Conf, CoC, U			No	Yes
BNZX	Benzodiazepines Conf, CoC, U			No	Yes

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
OXYCX	36300	Oxycodone Immunoassay Screen	Alphanumeric	ng/mL	19642-8
OXYCX	61728	Oxycodone-by LC-MS/MS	Alphanumeric	ng/mL	16249-5
OXYCX	42012	Noroxycodone-by LC-MS/MS	Alphanumeric	ng/mL	61425-5
OXYCX	35964	Oxymorphone-by LC-MS/MS	Alphanumeric	ng/mL	17395-5
OXYCX	42013	Noroxymorphone-by LC-MS/MS	Alphanumeric	ng/mL	90894-7
OXYCX	36023	Oxycodone Interpretation	Alphanumeric		69050-3
OXYCX	36028	Chain of Custody	Alphanumeric		77202-0
ADLTX	36121	Creatinine, U	Numeric	mg/dL	2161-8
ADLTX	36122	Specific Gravity	Numeric		5810-7
ADLTX	36123	pH	Numeric		2756-5
ADLTX	36124	Oxidants	Alphanumeric		58714-7
ADLTX	36125	Nitrites	Alphanumeric		32710-6
ADLTX	36126	Comment	Alphanumeric		48767-8

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
AMPHX	6538	Amphetamines Immunoassay Screen	Alphanumeric		19261-7
AMPHX	36128	Amphetamine-by LC-MS/MS	Alphanumeric	ng/mL	20410-7
AMPHX	36129	Phentermine-by LC-MS/MS	Alphanumeric	ng/mL	20557-5
AMPHX	36130	Methamphetamine-by LC-MS/MS	Alphanumeric	ng/mL	16235-4
AMPHX	36131	Pseudoephedrine/Ephedrine-by LC-MS/MS	Alphanumeric	ng/mL	58707-1
AMPHX	36132	MDA (Ecstasy metabolite)-by LC-MS/MS	Alphanumeric	ng/mL	20545-0
AMPHX	36133	MDMA (Ecstasy)-by LC-MS/MS	Alphanumeric	ng/mL	18358-2
AMPHX	36134	Amphetamines Interpretation	Alphanumeric		69050-3
AMPHX	36135	Chain of Custody	Alphanumeric		77202-0
COKEX	2903	Cocaine Immunoassay Screen	Alphanumeric		42241-0
COKEX	36162	Cocaine-by GC/MS	Alphanumeric	ng/mL	20519-5
COKEX	36163	Benzoyllecgonine-by GC/MS	Alphanumeric	ng/mL	16226-3
COKEX	36164	Cocaine Interpretation	Alphanumeric		69050-3
COKEX	36165	Chain of Custody	Alphanumeric		77202-0
BARBX	6668	Barbiturates Immunoassay Screen	Alphanumeric		70155-7
BARBX	36143	Butalbital-by GC/MS	Alphanumeric	ng/mL	16237-0
BARBX	36144	Amobarbital-by GC/MS	Alphanumeric	ng/mL	16239-6
BARBX	36145	Pentobarbital-by GC/MS	Alphanumeric	ng/mL	16240-4
BARBX	36146	Secobarbital-by GC/MS	Alphanumeric	ng/mL	16238-8
BARBX	36147	Phenobarbital-by GC/MS	Alphanumeric	ng/mL	16241-2
BARBX	36148	Barbiturates Interpretation	Alphanumeric		69050-3
BARBX	36149	Chain of Custody	Alphanumeric		77202-0
MTDNX	21105	Methadone Immunoassay Screen	Alphanumeric		70149-0
MTDNX	36208	EDDP-by GC-MS	Alphanumeric	ng/mL	58429-2
MTDNX	36209	Methadone-by GC-MS	Alphanumeric	ng/mL	16246-1
MTDNX	36210	Methadone Interpretation	Alphanumeric		69050-3
MTDNX	36211	Chain of Custody	Alphanumeric		77202-0
PCPX	6672	Phencyclidine Immunoassay Screen	Alphanumeric		19659-2
PCPX	36229	Phencyclidine-by GC/MS	Alphanumeric	ng/mL	16254-5
PCPX	36230	Phencyclidine Interpretation	Alphanumeric		69050-3
PCPX	36231	Chain of Custody	Alphanumeric		77202-0

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
THCX	2449	Carboxy-THC Immunoassay Screen	Alphanumeric		19415-9
THCX	618771	Delta-8 Carboxy-Tetrahydrocannabinol by LC-MS/MS	Alphanumeric	ng/mL	20521-1
THCX	36238	Delta-9 Carboxy-Tetrahydrocannabinol by LC-MS/MS	Alphanumeric	ng/mL	20521-1
THCX	36239	Carboxy-THC Interpretation	Alphanumeric		69050-3
THCX	36240	Chain of Custody	Alphanumeric		77202-0
OPATX	6541	Opiates Immunoassay Screen	Alphanumeric	ng/mL	70151-6
OPATX	36213	Codeine-by LC-MS/MS	Alphanumeric	ng/mL	16250-3
OPATX	42005	Dihydrocodeine-by LC-MS/MS	Alphanumeric	ng/mL	19448-0
OPATX	36214	Hydrocodone-by LC-MS/MS	Alphanumeric	ng/mL	16252-9
OPATX	42006	Norhydrocodone-by LC-MS/MS	Alphanumeric	ng/mL	61422-2
OPATX	36215	Hydromorphone-by LC-MS/MS	Alphanumeric	ng/mL	16998-7
OPATX	36216	Oxycodone-by LC-MS/MS	Alphanumeric	ng/mL	16249-5
OPATX	42007	Noroxycodone-by LC-MS/MS	Alphanumeric	ng/mL	61425-5
OPATX	36212	Oxymorphone-by LC-MS/MS	Alphanumeric	ng/mL	17395-5
OPATX	42008	Noroxymorphone-by LC-MS/MS	Alphanumeric	ng/mL	90894-7
OPATX	42009	Naloxone-by LC-MS/MS	Alphanumeric	ng/mL	77207-9
OPATX	36217	Morphine-by LC-MS/MS	Alphanumeric	ng/mL	16251-1
OPATX	36218	Opiates Interpretation	Alphanumeric		69050-3
OPATX	36219	Chain of Custody	Alphanumeric		77202-0
BNZX	608452	Benzodiazepines Immunoassay Screen	Alphanumeric		14316-4
BNZX	608280	Alprazolam by LC-MS/MS	Alphanumeric	ng/mL	59615-5
BNZX	608281	Alpha-Hydroxyalprazolam by LC-MS/MS	Alphanumeric	ng/mL	16348-5
BNZX	608282	Chlordiazepoxide by LC-MS/MS	Alphanumeric	ng/mL	20522-9
BNZX	608283	Clonazepam by LC-MS/MS	Alphanumeric	ng/mL	16229-7
BNZX	608284	7-aminoclonazepam by LC-MS/MS	Alphanumeric	ng/mL	51776-3
BNZX	608285	Diazepam by LC-MS/MS	Alphanumeric	ng/mL	16227-1
BNZX	608286	Nordiazepam by LC-MS/MS	Alphanumeric	ng/mL	16228-9
BNZX	608287	Midazolam by LC-MS/MS	Alphanumeric	ng/mL	16233-9
BNZX	608288	Alpha-Hydroxy Midazolam by LC-MS/MS	Alphanumeric	ng/mL	59590-0

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
BNZX	608289	Oxazepam by LC-MS/MS	Alphanumeric	ng/mL	16201-6
BNZX	608290	Temazepam by LC-MS/MS	Alphanumeric	ng/mL	20559-1
BNZX	608291	Clobazam by LC-MS/MS	Alphanumeric	ng/mL	59635-3
BNZX	608292	N-Desmethyclobazam by LC-MS/MS	Alphanumeric	ng/mL	97160-6
BNZX	608293	Flunitrazepam by LC-MS/MS	Alphanumeric	ng/mL	20528-6
BNZX	608294	7-aminoflunitrazepam by LC-MS/MS	Alphanumeric	ng/mL	51777-1
BNZX	608295	Flurazepam by LC-MS/MS	Alphanumeric	ng/mL	16231-3
BNZX	608296	2-Hydroxy Ethyl Flurazepam by LC-MS/MS	Alphanumeric	ng/mL	97159-8
BNZX	608297	Lorazepam by LC-MS/MS	Alphanumeric	ng/mL	17088-6
BNZX	608298	Prazepam by LC-MS/MS	Alphanumeric	ng/mL	17479-7
BNZX	608299	Triazolam by LC-MS/MS	Alphanumeric	ng/mL	16232-1
BNZX	608300	Alpha-Hydroxy Triazolam by LC-MS/MS	Alphanumeric	ng/mL	49876-6
BNZX	608301	Zolpidem by LC-MS/MS	Alphanumeric	ng/mL	72770-1
BNZX	608302	Zolpidem Phenyl-4-Carboxylic acid by LC-MS/MS	Alphanumeric	ng/mL	72768-5
BNZX	608449	Benzodiazepines Interpretation	Alphanumeric		69050-3
BNZX	608450	Chain of Custody	Alphanumeric		77202-0

Reference Values:

Negative

Screening cutoff concentrations

Amphetamines: 500 ng/mL

Barbiturates: 200 ng/mL

Benzodiazepines: 100 ng/mL

Cocaine (benzoylecgonine-cocaine metabolite): 150 ng/mL

Methadone metabolite: 300 ng/mL

Opiates: 300 ng/mL

Phencyclidine: 25 ng/mL

Tetrahydrocannabinol carboxylic acid: 50 ng/mL

This report is intended for use in clinical monitoring or management of patients. It is not intended for use in employment-related testing.