

### **Test Definition: FENR**

Fentanyl Screen with Reflex, Random, Urine

Reporting Title: Fentanyl Screen w/Reflex, U

Performing Location: Rochester

### **Ordering Guidance:**

For situations where chain of custody is required, a Chain of Custody Kit (T282) is available. For chain-of-custody testing, order FENTX / Fentanyl with Metabolite Confirmation, Chain of Custody, Urine.

For monitoring therapeutic drug levels, order FENTS / Fentanyl, Serum.

#### **Additional Testing Requirements:**

If urine creatinine is required or adulteration of the sample is suspected, the ADULT / Adulterants Survey, Random, Urine test should also be ordered.

### **Specimen Requirements:**

Supplies: Sarstedt Aliquot Tube, 5 mL (T914)

Collection Container/Tube: Clean, plastic urine collection container

Submission Container/Tube: Plastic, 5-mL tube

Specimen Volume: 5 mL Collection Instructions:

1. Collect a random urine specimen.

2. No preservative.

#### **Specimen Minimum Volume:**

2.5 mL

#### Forms:

If not ordering electronically, complete, print, and send a Therapeutics Test Request (T831) with the specimen.

| Specimen Type | Temperature              | Time     | Special Container |
|---------------|--------------------------|----------|-------------------|
| Urine         | Refrigerated (preferred) | 14 days  |                   |
|               | Frozen                   | 14 days  |                   |
|               | Ambient                  | 72 hours |                   |



## **Test Definition: FENR**

Fentanyl Screen with Reflex, Random, Urine

### **Result Codes:**

| Result ID | Reporting Name           | Туре         | Unit  | LOINC®  |
|-----------|--------------------------|--------------|-------|---------|
| 63060     | Fentanyl Screen, U       | Alphanumeric | ng/mL | 59673-4 |
|           | Also used by tests: FENS |              |       |         |

LOINC and CPT codes are provided by the performing laboratory.

### **Supplemental Report:**

No

### **Components:**

| Test ID | Reporting Name     | CPT Units | CPT Code | Always<br>Performed | Orderable<br>Separately |
|---------|--------------------|-----------|----------|---------------------|-------------------------|
| FENS    | Fentanyl Screen, U |           |          | Yes                 | Yes                     |

#### **CPT Code Information:**

80307

### **Reflex Tests:**

| Test ID | Reporting Name                | CPT Units | CPT Code | Always<br>Performed | Orderable<br>Separately |
|---------|-------------------------------|-----------|----------|---------------------|-------------------------|
| FENTU   | Fentanyl w/metabolite Conf, U |           |          | No                  | Yes                     |

### **Result Codes for Reflex Tests:**

| Test ID | Result ID | Reporting Name          | Туре         | Unit  | LOINC®  |
|---------|-----------|-------------------------|--------------|-------|---------|
| FENTU   | 45530     | Fentanyl by LC-MS/MS    | Alphanumeric | ng/mL | 58381-5 |
| FENTU   | 45531     | Norfentanyl by LC-MS/MS | Alphanumeric | ng/mL | 58383-1 |
| FENTU   | 45532     | Fentanyl Interpretation | Alphanumeric |       | 69050-3 |



# **Test Definition: FENR**

Fentanyl Screen with Reflex, Random, Urine

### **Reference Values:**

Negative

Screening cutoff concentration: 2 ng/mL