

**Test Definition: CLAV** 

Cortisol, Left Adrenal Vein, Serum

**Reporting Title:** Cortisol, LAV **Performing Location:** Rochester

## Specimen Requirements:

Collection Container/Tube: Red top

Specimen Volume: 1.5 mL

Submission Container/Tube: Plastic vial

Collection Instructions:

- 1. Morning (8 a.m.) and afternoon (4 p.m.) specimens are preferred.
- 2. Include time of collection
- 3. Centrifuge and aliquot serum into a plastic vial.

Additional Information: If multiple specimens are collected, send a separate order for each specimen.

#### **Specimen Minimum Volume:**

0.5 mL

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	28 days	

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
6346	Cortisol, LAV	Numeric	mcg/dL	2143-6

LOINC and CPT codes are provided by the performing laboratory.

## **Supplemental Report:**

No

#### **CPT Code Information:**

82533



# **Test Definition: CLAV**

Cortisol, Left Adrenal Vein, Serum

## **Reference Values:**

No established reference values