Wilms Tumor (WT-1) Immunostain, Technical Component Only

**Test Definition: WT1I** 

# **Reporting Title:** Wilm's Tumor (WT-1) IHC, Tech Only **Performing Location:** Rochester

MAYO CLINIC LABORATORIES

### Ordering Guidance:

This test includes only technical performance of the stain (no pathologist interpretation is performed). If diagnostic consultation by a pathologist is required order PATHC / Pathology Consultation.

#### Shipping Instructions:

Attach the green "Attention Pathology" address label (T498) and the pink Immunostain Technical Only label included in the kit to the outside of the transport container.

#### **Specimen Requirements:**

Supplies: Immunostain Technical Only Envelope (T693) Specimen Type: Tissue Container/Tube: Immunostains Technical Only Envelope Preferred: 2 Unstained positively charged glass slide (25- x 75- x 1-mm) per test ordered; sections 4-microns thick. Acceptable: Formalin-fixed, paraffin-embedded tissue block

#### Forms:

If not ordering electronically, complete, print, and send a Immunohistochemical (IHC)/In Situ Hybridization (ISH) Stains Request (T763) with the specimen.

Specimen Type	Temperature	Time	Special Container
TECHONLY	Ambient (preferred)		
	Refrigerated		

# **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
70894	Wilm's Tumor (WT-1) IHC, Tech Only	Alphanumeric		Bill only; no re

LOINC and CPT codes are provided by the performing laboratory.

# Supplemental Report:

No



# **Reflex Tests:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
ΙΗΤΟΑ	A IHC Additional, Tech Only			No	No
IHTOI IHC Initial, Tech Only				No	No