

Test Definition: SUBIF

Immunoglobulin G (IgG) Subtypes Immunofluorescence, Tissue

Reporting Title: IGG Subtypes Immunofluorescence

Performing Location: Rochester

Ordering Guidance:

If additional interpretation/analysis is needed, request PATHC / Pathology Consultation along with this test and send the corresponding renal pathology light microscopy and immunofluorescence (IF) slides (or IF images on a CD), electron microscopy images (prints or CD), and the pathology report.

Shipping Instructions:

- 1. Advise shipping specimens in Styrofoam transportation coolers to avoid extreme hot or cold temperatures to ensure specimens are received at required specimen stability temperature.
- 2. Attach the green pathology address label included in the kit to the outside of the transport container.

Necessary Information:

A pathology/diagnostic report is required.

Specimen Requirements:

Preferred: Frozen tissue

Supplies: Renal Biopsy Kit (T231) Specimen Type: Kidney tissue

Container/Tube: Renal Biopsy Kit, Zeus/Michel's, Frozen

Specimen Volume: Entire specimen

Collection Instructions: Collect specimens according to the instructions in Renal Biopsy Procedure for Handling Tissue

for Light Microscopy, Immunofluorescent Histology, and Electron Microscopy.

Additional Information: If standard immunoglobulin and complement immunofluorescence has already been performed,

submit the residual frozen tissue (must contain glomeruli) on dry ice.

Acceptable: Frozen tissue

Slides: 4 frozen tissue unstained positively-charged glass slides (25- x 75- x 1-mm) per test ordered; sections 4-microns

thick.

Forms:

- 1. Renal Biopsy Patient Information
- 2. If not ordering electronically, complete, print, and send a Renal Diagnostics Test Request (T830) with the specimen.

Specimen Type	Temperature	Time	Special Container
Special	Frozen (preferred)		



Test Definition: SUBIF

Immunoglobulin G (IgG) Subtypes Immunofluorescence, Tissue

Ambient	
Refrigerated	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
BA0271	Interpretation	Alphanumeric		50595-8
BA0272	Participated in the Interpretation	Alphanumeric		No LOINC Needed
BA0275	Report electronically signed by	Alphanumeric		19139-5
BA0276	Addendum	Alphanumeric		35265-8
BA0274	Gross Description	Alphanumeric		22634-0
BA0273	Material Received	Alphanumeric		22633-2
71618	Disclaimer	Alphanumeric		62364-5
71850	Case Number	Alphanumeric		80398-1

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

88346-primary IF 88350-if additional IF