

 MAYO CLINIC
 Test Definition: LBNA∠

 LABORATORIES
 Epstein Barr Nuclear Antigen 2 (EBNA2) Immunostain, Technical

Component Only

Reporting Title: EBNA2 IHC, Tech Only Performing Location: Rochester

Ordering Guidance:

This test includes only technical performance of the stain (no pathologist interpretation is performed). If diagnostic consultation by a pathologist is required order PATHC / Pathology Consultation.

Shipping Instructions:

Attach the green "Attention Pathology" address label (T498) and the pink Immunostain Technical Only label included in the kit to the outside of the transport container.

Specimen Requirements:

Supplies: Immunostain Technical Only Envelope (T693)

Specimen Type: Tissue

Container/Tube: Immunostain Technical Only Envelope

Preferred: 2 Unstained positively charged glass slide (25- x 75- x 1-mm) per test ordered; sections 4-microns thick Acceptable: Formalin-fixed, paraffin-embedded tissue block

Forms:

If not ordering electronically, complete, print, and send an Immunohistochemical (IHC)/In Situ Hybridization (ISH) Stains Request (T763) with the specimen.

Specimen Type	Temperature	Time	Special Container
TECHONLY	Ambient (preferred)		
	Refrigerated		

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
71494	EBNA2 IHC, Tech Only	Alphanumeric		Bill only; no re

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:



 MAYOCLINIC
 Test Definition: EBNA2

 LABORATORIES
 Epstein Barr Nuclear Antigen 2 (EBNA2) Immunostain, Technical

 Component Only
 Component Only

88342-TC, primary 88341-TC, if additional IHC

Reflex Tests:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
IHTOI	IHC Initial, Tech Only			No	No
IHTOA	OA IHC Additional, Tech Only			No	No