Meconium Methadone Screen with Reflex Confirmation

Test Definition: FFMTH

Reporting Title: Meconium Methadone Scr/Confirm **Performing Location:** Medtox Laboratories, Inc.

Specimen Requirements:

Container/Tube: Sterile collection container

MAYO CLINIC LABORATORIES

Specimen Volume: 3 g

Collection Instructions: Collect specimen into the same sterile collection container until 3 g of meconium have been collected or until the first milk stool appears. When at least 3 g of meconium has been collected, tightly screw on the cap of the collection vial and send specimen ambient.

Note: Specimens from different voids may be pooled if necessary.

Specimen Minimum Volume:

1 g

Specimen Type	Temperature	Time	Special Container
Meconium	Ambient (preferred)	14 days	
	Frozen	365 days	
	Refrigerated	14 days	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
Z5603	Methadone	Alphanumeric		43595-8

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

80307



Reflex Tests:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
FMTHM	Meconium Methadone Confirm			No	No

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Туре	Unit	LOINC®
FMTHM	Z5604	Methadone	Alphanumeric		29159-1
FMTHM	Z5605	Methadone Metabolite (EDDP)	Alphanumeric		73562-1

Reference Values:

Initial presumptive testing by immunoassay at the following testing threshold: Methadone: 50 ng/gm

Presumptive positive confirmed by definitive chromatography mass spectrometry (GC/MC or LC/MS-MS)