

Reporting Title: Aspergillus Ab, Qn, DID**Performing Location:** LabCorp Burlington**Specimen Requirements:**

Specimen Type: Serum

Container/Tube: Red top or SST

Specimen Volume: 1 mL

Collection Instructions: Draw blood in a plain red-top tube, serum gel tube(s) is also acceptable. Separate serum immediately after coagulation (30 minutes) to prevent hemolysis. Send 1 mL of serum frozen in a plastic vial.

NOTE: Patient should be fasting for eight hours to avoid lipemic sample interference.

Specimen Minimum Volume:

0.5 mL

Specimen Type	Temperature	Time	Special Container
Serum	Frozen (preferred)	14 days	
	Refrigerated	72 hours	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z5651	Aspergillus fumigatus	Alphanumeric		23821-2
Z5652	Aspergillus flavus	Alphanumeric		23820-4
Z5653	Aspergillus niger	Alphanumeric		10894-4

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

86606 x 3

Reference Values:

Aspergillus fumigatus	Neg: <1:1
Aspergillus flavus	Neg: <1:1
Aspergillus niger	Neg: <1:1