

**Test Definition: FVENE** 

Venison IgE

Reporting Title: Venison IgE

Performing Location: Eurofins Viracor

# **Specimen Requirements:**

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

# **Specimen Minimum Volume:**

0.5 mL

Specimen Type	Temperature	Time	Special Container	
Serum	Refrigerated (preferred)	28 days		
	Frozen	365 days		
	Ambient	28 days		

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
Z5660	Venison IgE	Alphanumeric		7759-4
Z5661	Class	Alphanumeric		102840-6

LOINC and CPT codes are provided by the performing laboratory.

### **Supplemental Report:**

No

# **CPT Code Information:**

86003



**Test Definition: FVENE** 

Venison IgE

### **Reference Values:**

<0.35 kU/L