

Reporting Title: ANA, Synovial Fluid**Performing Location:** Esoterix Endocrinology**Specimen Requirements:**

Specimen Type: Synovial fluid

Collection Container/Tube: Sterile container

Submission Container/Tube: Plastic vial

Specimen Volume: 1 mL

Collection Instructions:

1. Collect synovial fluid and transfer into a plastic vial and ship refrigerate.
2. Indicate the specimen source.

Specimen Minimum Volume:

0.3 mL

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	14 days	
	Frozen	60 days	
	Ambient	7 days	

Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Type	Reportable
FFANA	Z5710	ANA Specimen Type	Plain Text	No

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z5711	ANA, Synovial Fluid	Alphanumeric		In Process

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

86038

86039 (if appropriate)

Reflex Tests:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
FPANA	ANA Titer and Pattern Synovial Fluid			No	No

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
FPANA	Z5729	Homogeneous Pattern	Alphanumeric		13067-4
FPANA	Z5730	Nucleolar Pattern	Alphanumeric		13067-4
FPANA	Z5731	Speckled Pattern	Alphanumeric		13067-4
FPANA	Z5732	Centromere Pattern	Alphanumeric		13067-4
FPANA	Z5733	Spindle Apparatus Pattern	Alphanumeric		13067-4
FPANA	Z5734	Nuclear Membrane Pattern	Alphanumeric		13067-4
FPANA	Z5735	Midbody Pattern	Alphanumeric		13067-4
FPANA	Z5736	Nuclear Dot Pattern	Alphanumeric		13067-4
FPANA	Z5737	PCNA Pattern	Alphanumeric		13067-4
FPANA	Z5738	Centriole Pattern	Alphanumeric		13067-4
FPANA	Z5712	Note:	Alphanumeric		13067-4

Reference Values:

ANA Titer: <1:10
ANA Pattern: No Pattern