

Reporting Title: SeroNeg RAdx3 Profile**Performing Location:** Esoterix Endocrinology**Specimen Requirements:**

Collection Container/Tube: 10 mL Red

Submission Container/Tube: Plastic vial

Specimen Volume: 3 mL

Acceptable: SST

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube is acceptable. Spin down and send 3 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume:

3 mL (volume does NOT allow for repeat testing)

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	14 days	
	Frozen	60 days	
	Ambient	7 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z5779	14.3.3 ETA, Rheum. Arthritis	Alphanumeric		75880-5
Z5780	Anti-CEP-1 Ab, IgG	Alphanumeric		96561-6
Z5781	Anti-Sa Ab, IgG	Alphanumeric		54022-9

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

83520

83516 x2

Reference Values:

Interpretation for:

14.3.3 ETA, Rheum. Arthritis:

Reference Range: <0.20 ng/mL

Negative: <0.20 ng/mL

Weak Positive: 0.20-0.39 ng/mL

Moderate Positive: 0.40-0.79 ng/mL

Strong Positive: >0.79 ng/mL

Interpretation for:

Anti-CEP-1 Ab, IgG and Anti-Sa Ab, IgG:

Reference Range: <20

Negative: <20 units

Weak Positive: 20-39 units

Moderate Positive: 40-80 units

Strong Positive: >80 units