

Reporting Title: Anti-MDA-5 Ab (CADM-140)**Performing Location:** Esoterix Endocrinology**Ordering Guidance:**

This test is a component of the Myomarker 3 Profile and the Myomarker 3 Profile Plus.

Specimen Requirements:

Collection Container/Tube: Red

Submission Container/Tube: Plastic vial

Specimen Volume: 1 mL

Acceptable: SST

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube is acceptable. Spin down and send 1 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume:

0.3 mL (volume does NOT allow for repeat testing)

| Specimen Type | Temperature | Time | Special Container |
|---------------|--------------------------|---------|-------------------|
| Serum | Refrigerated (preferred) | 14 days | |
| | Frozen | 60 days | |
| | Ambient | 7 days | |

Result Codes:

| Result ID | Reporting Name | Type | Unit | LOINC® |
|-----------|--------------------------|--------------|------|--------------|
| FAMDA | Anti-MDA-5 Ab (CADM-140) | Alphanumeric | | Not Provided |

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

83520

Reference Values:

Reference Range: <20

Interpretation:

Negative: <20 units

Weak Positive: 20-39 units

Moderate Positive: 40-80 units

Strong Positive: >80 units