

**Reporting Title:** Anti-SAE1 Ab, IgG**Performing Location:** Esoterix Endocrinology**Specimen Requirements:**

Collection Container/Tube: 5 mL Red

Submission Container/Tube: Plastic vial

Specimen Volume: 1 mL

Acceptable: SST

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube is acceptable. Spin down and send 1 mL of serum refrigerated in a plastic vial.

**Specimen Minimum Volume:**

0.5 mL (volume does NOT allow for repeat testing)

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	14 days	
	Frozen	60 days	
	Ambient	7 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
FAS1A	Anti-SAE1 Ab, IgG	Alphanumeric		101220-2

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

83516

---

**Reference Values:**

Reference Range: <20

Interpretation:

Negative: <20 units

Weak Positive: 20-39 units

Moderate Positive: 40-80 units

Strong Positive: >80 units