

Test Definition: FANIC

Anti-Nuclear Antibodies by Indirect Fluorescent Antibody (IFA), Cerebrospinal Fluid

Reporting Title: ANA by IFA, CSF

Performing Location: Esoterix Endocrinology

Specimen Requirements:

Specimen Type: Cerebrospinal fluid (CSF) Collection Container/Tube: Sterile container Submission Container/Tube: Plastic vial

Specimen Volume: 2 mL

Collection Instructions: Collect cerebrospinal fluid (CSF) in a sterile container and ship refrigerate.

Specimen Minimum Volume:

0.5 mL (Note: this volume does NOT allow for repeat testing)

Specimen Type	Temperature	Time	Special Container
CSF	Refrigerated (preferred)	14 days	
	Frozen	60 days	
	Ambient	7 days	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
FANIC	ANA by IFA, CSF	Alphanumeric		63294-3

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

86038

86039 (if appropriate)



Test Definition: FANIC

Anti-Nuclear Antibodies by Indirect Fluorescent Antibody (IFA), Cerebrospinal Fluid

Reflex Tests:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
FANTC	ANA Titer and Pattern CSF			No	No

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Туре	Unit	LOINC®
FANTC	Z5826	Homogeneous Pattern	Alphanumeric		13067-4
FANTC	Z5827	Nucleolar Pattern	Alphanumeric		13067-4
FANTC	Z5828	Speckled Pattern	Alphanumeric		13067-4
FANTC	Z5829	Centromere Pattern	Alphanumeric		13067-4
FANTC	Z5830	Spindle Apparatus Pattern	Alphanumeric		13067-4
FANTC	Z5831	Nuclear Membrane Pattern	Alphanumeric		13067-4
FANTC	Z5832	Midbody Pattern	Alphanumeric		13067-4
FANTC	Z5833	Nuclear Dot Pattern	Alphanumeric		13067-4
FANTC	Z5834	PCNA Pattern	Alphanumeric		13067-4
FANTC	Z5835	Centriole Pattern	Alphanumeric		13067-4
FANTC	Z5836	Note:	Alphanumeric		Not Provided

Reference Values:

ANA by IFA, CSF: Negative ANA Pattern: No Pattern