

**Reporting Title:** Thrombin Antithrombin Complex**Performing Location:** Esoterix Coagulation**Specimen Requirements:**

Specimen Type: Citrated plasma

Collection Container/Tube: Light-blue top (citrate)

Specimen Volume: 2 mL

Collection Instructions: Draw blood in a light blue-top (Sodium citrate) tube(s). Spin down and send 2 mL citrated plasma frozen in a plastic vial.

Note: Separate specimens must be submitted when multiple tests are ordered.

**Specimen Minimum Volume:**

1 mL

Specimen Type	Temperature	Time	Special Container
Plasma Na Cit	Frozen	180 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
FFTAC	Thrombin Antithrombin Complex	Alphanumeric		14182-0

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

83520

---

**Reference Values:**

<4.3 ng/mL

Pre-analytical conditions such as a difficult draw may spuriously increase test results.