

Reporting Title: Aspergillus IgG Precipitins Panel**Performing Location:** Eurofins Viracor**Specimen Requirements:**

Specimen Type: Serum

Container/Tube: Red or SST

Specimen Volume: 1 mL

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube is acceptable. Spin down and send 1 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume:

0.5 mL

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	28 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z5917	Aspergillus amstelodami/glaucus	Alphanumeric		Not Provided
Z5918	Aspergillus flavus	Alphanumeric		Not Provided
Z5919	Aspergillus fumigatus Mix	Alphanumeric		Not Provided
Z5920	Aspergillus nidulans	Alphanumeric		Not Provided
Z5921	Aspergillus niger	Alphanumeric		Not Provided
Z5922	Aspergillus versicolor	Alphanumeric		Not Provided

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

86331 x 6

Reference Values:

Negative