

# **Test Definition: FIXCF**

Immunofixation, CSF

Reporting Title: Immunofixation, CSF

Performing Location: Quest Diagnostics Nichols Institute

### **Specimen Requirements:**

Specimen Type: CSF

Container/Tube: sterile screw cap container

Specimen Volume: 6 mL

Collection Instructions: Collect 6 mL of spinal fluid (CSF) in a sterile screw cap container. Ship frozen.

## **Specimen Minimum Volume:**

4.5 mL

Specimen Type	Temperature	Time	Special Container
CSF	Frozen (preferred)	30 days	
	Refrigerated	10 days	
	Ambient	6 days	

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
FIXCF	Immunofixation, CSF	Alphanumeric		13174-8

LOINC and CPT codes are provided by the performing laboratory.

### **Supplemental Report:**

Referral

#### **CPT Code Information:**

86335



# **Test Definition: FIXCF**

Immunofixation, CSF

### **Reference Values:**

No abnormal bands are present on immunofixation.