

**Reporting Title:** CIC, C3 fragments  
**Performing Location:** ARUP Laboratories

**Specimen Requirements:**

Specimen Type: Serum

Container/Tube: Red Top/SST acceptable

Specimen Volume: 1 mL

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube(s), is acceptable. Allow complete clotting of red blood cells (up to 1 hour), then separate from cells within 30 minutes and freeze immediately. Send 1 mL of serum frozen in a plastic vial.

Note: Critical frozen. Separate specimens must be submitted when multiple tests are ordered.

**Specimen Minimum Volume:**

0.5 mL

Specimen Type	Temperature	Time	Special Container
Serum	Frozen	30 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
Z6125	CIC, C3 fragments	Alphanumeric		5228-2

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

86332

---

**Reference Values:**

Less than or equal to 15 ug Eq/mL