

Reporting Title: VIP, Plasma**Performing Location:** Quest Diagnostics Nichols Institute**Specimen Requirements:**

Specimen Type: Plasma

Collection Container/Tube: Lavender top (EDTA)

Submission Container/Tube: Preservative-free plastic vial

Specimen Volume: 1 mL

Collection Instructions: Draw blood in an EDTA (lavender top) tube(s). Spin down and freeze immediately. Send 1 mL of plasma frozen in a preservative-free plastic vial.

Note: Fasting is preferred.

Specimen Minimum Volume:

0.75 mL

Specimen Type	Temperature	Time	Special Container
Plasma EDTA	Frozen	180 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
FVIPL	VIP, Plasma	Alphanumeric		3125-2

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

84586

Reference Values:

<78 pg/mL